

Multi-Use, Ultrasound Table with Stirrups



MODEL 80069

- Three-piece top with pneumatic adjustable backrest, manual footrest, and stirrups
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Trendelenburg/Reverse Trendelenburg plus or minus 15 °
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Industry leading 5-year, limited warranty
- Optional side rails & casters available
- 34" wide model available
- Model 80069: 500 lbs load capacity under normal use
- Model 80069-X: 600 lbs load capacity under normal use



Trendelenburg



Reverse Trendelenburg




Optional Side Rails & Casters

80069 SPECIFICATIONS

	length	width	height range
overall	76"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
backrest	38"	27"	
footrest	12"	18 ¹ / ₂ "	

80069-X Extra-Wide 34" width model

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



Upholstery Colors



OPTIONS

Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set

080	length	width
	92"	2"
	233.67 cm	5 cm



Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

098



Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087



Paper Cutter

- Clear plastic with elastic strap and buckle

040

