

Power, Imaging Table with Three-Section Top and Drop Window



MODEL 80073

- Three-piece top with pneumatic adjustable backrest, manual footrest and window drop
- Footrest is adjustable from either side and raises up to 30° from flat position
- Optional side rails & casters (see options on page 2)
- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Industry leading 5-year, limited warranty
- 34" wide model available
- Model 80073: 500 lbs load capacity under normal use
- Model 80073-X: 600 lbs load capacity under normal use

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

80073 SPECIFICATIONS

	length	width	height range
overall	76"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
backrest	38"	27"	
footrest	12"	27"	
window	10"	10"	

80073-X Extra-Wide 34" width model



Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

OPTIONS

Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set

080	length	width
	92"	2"
	233.67 cm	5 cm



Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

098



Power Backrest

- Smooth, power control raises to 80° and lowers to flat

011



Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087



Paper Cutter

- Clear plastic with elastic strap and buckle

040

