

Panel Leg Series, Treatment Table with Shelving



MODELS: 91030-27 & 91030-30


- All laminate construction
- Sleek panel leg design
- Laminate shelf and tiered storage compartment
- Adjustable backrest with 6 settings of approximately: 27°, 36°, 43°, 50°, 55° & 60°
- Scratch resistant topcoat finish
- Ships flat for easy handling and extra savings
- Fast and easy assembly
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 500 lbs. load capacity under normal use (181.43 kg)

SPECIFICATIONS

91030-27	72" L x 27" W x 31" H
91030-30	72" L x 30" W x 31" H



Prop 65 Warning




WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



034– Paper Dispenser & Cutter Combo



040– PaperCutter



030– Paper Dispenser



088–Flat Top



009–Nose Hole

Laminate Colors



Gray



Maple



Dark Cherry

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black