

Portrait™ VSM Vital Signs Monitor supplies and accessories

5975002-EN

Revision 5

© 2023 GE HealthCare.

All specifications are subject to change without notice.

This document is for software version Portrait VSM Software V1.0.

About this list

This list indicates the supplies and accessories that are approved, specified and available for use with the patient monitoring system.

To ensure patient safety, use only supplies and accessories listed in this document or recommended by GE HealthCare.

Refer to the user manual for complete safety information.

Please note that some products are not available worldwide. You can check the availability with your local GE HealthCare representative.

NOTE

Some of the recommended accessories may have to be ordered from the OEM manufacturer.

Naming conventions

For technical documentation purposes this list refers to the legal manufacturers as follows:

- Covidien LLC.: Nellcor
- GE Medical Systems *Information Technologies*, Inc.: GEMS-IT
- GE Healthcare Finland Oy.: GE Finland
- Masimo Corporation: Masimo
- Graphic Controls corporation, now a subsidiary of Nissha Printing Co., Ltd.: GRAPHIC CONTROLS
- Welch Allyn, Inc.: Welch Allyn
- Exergen Corporation: Exergen
- HeTaida Technology company: HeTaida

Ordering information

Contact your local GE HealthCare representative to order supplies or accessories.

Accessories

SpO₂ accessories

GE TruSignal SpO₂ accessories

Table 1 Cable, GE Finland

GE HealthCare Part Number	Manufacturer Part Number	Description
TS-G3	TS-G3	TruSignal interconnect cable with GE connector, 3 m/10 ft

Table 2 Reusable sensors, Non-sterile, GE Finland

GE HealthCare Part Number	Manufacturer Part Number	Description
TS-F-D	TS-F-D	TruSignal Finger sensor, adult-pediatric, 1 m/3.3 ft, ≥20 kg (≥44 lbs)
TS-E-D	TS-E-D	TruSignal Ear Sensor, adult-pediatric, 1 m/3.3 ft, ≥10 kg (≥22 lbs)
TS-W-D	TS-W-D	TruSignal Wrap sensor, adult-pediatric, 1 m/3.3 ft, ≥3 kg (≥6 lbs)
TS-SE-3	TS-SE-3	TruSignal Sensitive Skin sensor, adult-pediatric-infant-neonate, 1 m/3.3 ft, 3 pcs
TS-F2-GE	TS-F2-GE	TruSignal Integrated Finger sensors with GE connector, 2 m/6.6 ft
TS-F4-GE	TS-F4-GE	TruSignal Integrated Finger sensors with GE connector, 4 m/13 ft
TS-E2-GE	TS-E2-GE	TruSignal Integrated Ear sensors with GE connector, 2 m/6.6 ft
TS-E4-GE	TS-E4-GE	TruSignal Integrated Ear sensors with GE connector, 4 m/13 ft
TS-SA4-GE	TS-SA4-GE	TruSignal Integrated FingerTip sensor with GE connector, 4 m/13 ft, >30 kg (>66 lbs)
TS-SA-D	TS-SA-D	TruSignal FingerTip sensor with D connector, 1 m/3.3 ft, >30 kg (>66 lbs)
TS-SP-D	TS-SP-D	TruSignal PediTip sensor with D connector, 1 m/3.3 ft, 15 to 30 kg (33 to 66 lbs)
TS-SP3-GE	TS-SP3-GE	TruSignal Intergrated PediTip sensor with GE connector, 3m/10 ft, 15 to 30 kg (33 to 66 lbs)

Table 3 Adhesive (single patient use) sensors, Non-sterile, GE Finland

GE HealthCare Part Number	Manufacturer Part Number	Description
TS-AF-10	TS-AF-10	TruSignal AllFit sensor, adult-pediatric-infant-neonate, 0.5 m/1.6 ft, 10 pcs
TS-AF-25	TS-AF-25	TruSignal AllFit sensor, adult-pediatric-infant-neonate, 0.5 m/1.6 ft, 25 pcs
TS-AP-10	TS-AP-10	TruSignal adult/pediatric sensors, 0.3 m/1 ft, 20 to 50 kg (44 to 110 lbs) or ≥50 kg (≥110 lbs), 10 pcs

Table 3 Adhesive (single patient use) sensors, Non-sterile, GE Finland (Table continued)

GE HealthCare Part Number	Manufacturer Part Number	Description
TS-AP-25	TS-AP-25	TruSignal adult/pediatric sensors, 0.3 m/1 ft, 20 to 50 kg (44 to 110 lbs) or ≥50 kg (≥110 lbs), 25 pcs
TS-AAW-10	TS-AAW-10	TruSignal Adult Adhesive Wrap Sensor, ≥20 kg (≥44 lbs), 10 pieces
TS-AAW-25	TS-AAW-25	TruSignal Adult Adhesive Wrap Sensor, ≥20 kg (≥44 lbs), 25 pieces
TS-PAW-10	TS-PAW-10	TruSignal Pediatric Adhesive Wrap Sensor, 3 to 20 kg (6.6 to 44 lbs), 10 pieces
TS-PAW-25	TS-PAW-25	TruSignal Pediatric Adhesive Wrap Sensor, 3 to 20 kg (6.6 to 44 lbs), 25 pieces

Masimo SpO₂ accessories

Table 4 Cables

GE HealthCare Part Number	Manufacturer Part Number	Masimo SpO ₂ Accessory Description
MASIMO-RD-4085	4085	Interconnect Cable for Masimo RD SET [®] , GE-12, 3.6 m/12 ft.
2027263-002	2016	Interconnect Cable for Masimo LNCS, LNC-10, 3 m/10 ft.

Table 5 Reusable sensors, Non-sterile, Masimo

GE HealthCare Part Number	Manufacturer Part Number	Masimo SpO ₂ Accessory Description
MASIMO-RD-4050	4050	Finger Sensor, Adult, RD SET DCI, > 30 kg (> 66 lbs)
MASIMO-RD-4051	4051	Finger Sensor, Pediatric, RD SET DCI-P, 10-50 kg (22-110 lbs)
2027258-001	1863	Finger Sensor, Adult, LNCS DC-I, > 30 kg (> 66 lbs)
2027259-001	1864	Finger Sensor, Pediatric, LNCS DC-IP, 10-50 kg (22-110 lbs)
2027261-001	1895	Ear Sensor, Adult, LNCS TC-I, > 30 kg (> 66 lbs)

Table 6 Adhesive (single patient use) sensors, Non-sterile, Masimo

GE HealthCare Part Number	Manufacturer Part Number	Masimo SpO ₂ Accessory Description
MASIMO-RD-4000	4000	Adhesive sensor, Adult, RD SET ADT, 20/box
MASIMO-RD-4001	4001	Adhesive sensor, Pediatric, RD SET PDT, 20/box
MASIMO-RD-4002	4002	Adhesive sensor, Infant, RD SET INF, 20/box
MASIMO-RD-4003	4003	Adhesive sensor, Neonatal, RD SET NEO, 20/box
MASIMO-RD-4004	4004	Adhesive sensor, Neonatal/pediatric, RD SET NEO PT, 20/box
2027253-001	1859	Finger Sensor, Adult, LNCS-Adtx, > 30 kg (> 66 lbs), 20/box
2027254-001	1860	Finger Sensor, Pediatric, LNCS-Pdtx, 10-50 kg (22-110 lbs), 20/box
2051068-001	2328	Adhesive sensor, LNCS Inf, Infant (20/box)
2051069-001	2329	Adhesive sensor, LNCS Neo, Neonatal/Adult (20/box)
2051070-001	2330	Adhesive sensor, LNCS Neo, Neonatal Preterm (20/box)

Nellcor accessories

Table 7 Cables, GE Finland

GE HealthCare Part Number	Manufacturer Part Number	Nellcor SpO ₂ Accessory Description
2021406-001	2021406-001	Interconnect Cable for Nellcor OxiMax 3 m/10 ft
2021406-002	2021406-002	Interconnect Cable for Nellcor OxiMax 1.2 m/4 ft

Table 8 Reusable sensors, Non-sterile, Nellcor

GE HealthCare Part Number	Manufacturer Part Number	Nellcor SpO ₂ Accessory Description
407705-006	DS-100A	Finger Sensor, Adult, DS100A, >40 kg (>88 lbs)
70124021	DS-100A	Finger Sensor, Adult, DS100A, >40 kg (>88 lbs)
414248-001	OXI-P/I	Toe Oxiband Sensor, Pediatric/Infant, OXI-P/I, 3-40 kg (6.6-88 lbs)

Table 9 Adhesive (single patient use) sensors, Sterile, Nellcor

GE HealthCare Part Number	Manufacturer Part Number	Nellcor SpO ₂ Accessory Description
70124027	MAXAI	Finger Sensor, Adult, MAX-A, >30 kg (>66 lbs), 24/box
70124022	MAXPI	Finger Sensor, Pediatric, MAX-P, 3-50 kg (6.6-110 lbs), 24/box
70124026	MAXII	Toe Sensor, Infant, MAX-I, 3-20 kg (6.6-44 lbs), 24/box
70124032	MAXNI	Toe/Foot Sensor, Neonatal/Adult, MAX-N, <3 kg or >40 kg (<6.6 or >88 lbs) , 24/box

Non-invasive blood pressure accessories

Table 10 Air hoses, GEMS-IT

GE HealthCare Part Number	Manufacturer Part Number	NIBP Accessory Description
2058203-001	2058203-001	Adult-Pediatric, Dinaclick Connector, 1.2 m/4 ft
2058203-002	2058203-002	Adult-Pediatric, Dinaclick Connector, 3.6 m/12 ft
2058204-002	2058204-002	Neo-Snap Connector, 2.4 m/8 ft
2058204-003	2058204-003	Neo-Snap Connector, 3.6 m/12 ft
2058203-002-A	2058203-002-A	Adult-Pediatric, Dinaclick Connector, 3.6M/12ft Specific Labeling
2058204-003-A	2058204-003-A	Neo-Snap Connector, 3.6M/12ft Specific Labeling

Table 11 Assortment pack, blood pressure cuffs, GEMS-IT

GE HealthCare Part Number	Manufacturer Part Number	NIBP Accessory Description
2059301-501	2059301-501	ASSORTMENT PACK DURA-CUF SMALL ADULT, ADULT, LARGE ADULT 2T 80369-5 DINAClick, 3/ PK
2059301-502	2059301-502	ASSORTMENT PACK DURA-CUF PEDIATRIC 2T DINAClick, 80369-5, 6/ PK

Table 11 Assortment pack, blood pressure cuffs, GEMS-IT (Table continued)

GE HealthCare Part Number	Manufacturer Part Number	NIBP Accessory Description
2059301-503	2059301-503	Assortment pack, DURA-CUF, Adult, 2T DinaClick, 6/PK
2059303-501	2059303-501	ASSORTMENT PACK SOFT-CUF ADULT 2T CLICK, 80369-5, 20/ PK
2059303-502	2059303-502	ASSORTMENT PACK SOFT-CUF PEDIATRIC 2T DINAClick, 80369-5, 20/ PK
2059304-501	2059304-501	ASSORTMENT PACK CLASSIC-CUF ADULT 2T CLICK, 80369-5, 20/ PK
2059304-502	2059304-502	ASSORTMENT PACK CLASSIC-CUF SMALL ADULT, ADULT, LARGE ADULT 2T DINAClick, 80369-5, 3/PK
2059305-002	2059305-002	Assortment pack; SOFT-CUF, NEONATAL 2 TB NEO SNAP, NEO1-2, NEO2-3, NEO3-5, NEO4-5, NEO5-5
2059306-002	2059306-002	Assortment pack; CLASSIC-CUF Neonatal; 2T SNAP; 20/PK
2059306-003	2059306-003	Assortment pack; CLASSIC-CUF NEO 3, NEO 4, NEO 5 2T NEOSNAP

Table 12 DURA-CUF® blood pressure cuffs, dual-tube, DinaClick connector, GEMS-IT

GE HealthCare Part Number	Manufacturer Part Number	NIBP Accessory Description
DUR-A1-2A-5S	DUR-A1-2A-5S	DURA-CUF, SMALL ADULT, DINAClick, 17 - 25 CM, 80369-5, 5/ BOX
DUR-A1-2A-L-5S	DUR-A1-2A-L-5S	DURA-CUF, SMALL ADULT LONG, DINAClick, 17 - 25 CM, 80369-5, 5/ BOX
DUR-A2-2A-5S	DUR-A2-2A-5S	DURA-CUF, ADULT, DINAClick, 23 - 33 CM, 80369-5, 5/ BOX
DUR-A2-2A-L-5S	DUR-A2-2A-L-5S	DURA-CUF, ADULT LONG, DINAClick, 23 - 33 CM, 80369-5, 5/ BOX
DUR-A3-2A-5S	DUR-A3-2A-5S	DURA-CUF, LARGE ADULT, DINAClick, 31 - 40 CM, 80369-5, 5/ BOX
DUR-A3-2A-L-5S	DUR-A3-2A-L-5S	DURA-CUF, LARGE ADULT LONG, DINAClick, 31 - 40 CM, 80369-5, 5/ BOX
DUR-P1-2A-5S	DUR-P1-2A-5S	DURA-CUF, INFANT, DINAClick, 08 - 13 CM, 80369-5, 5/ BOX
DUR-P2-2A-5S	DUR-P2-2A-5S	DURA-CUF, CHILD, DINAClick, 12 - 19 CM, 80369-5, 5/ BOX
DUR-P1-2A-A	DUR-P1-2A-A	DURA-CUF, INFANT, DINAClick, 08 - 13 CM, 80369-5, 5/ BOX Specific Labeling
DUR-P2-2A-A	DUR-P2-2A-A	DURA-CUF, CHILD, DINAClick, 12 - 19 CM, 80369-5, 5/ BOX Specific Labeling
DUR-A1-2A-A	DUR-A1-2A-A	DURA-CUF, SMALL ADULT, DINAClick, 17 - 25 CM, 80369-5, 5/ BOX Specific Labeling
DUR-A2-2A-A	DUR-A2-2A-A	DURA-CUF, ADULT, DINAClick, 23 - 33 CM, 80369-5, 5/ BOX Specific Labeling
DUR-A2-2A-L-A	DUR-A2-2A-L-A	DURA-CUF, ADULT LONG, DINAClick, 23 - 33 CM, 80369-5, 5/ BOX Specific Labeling
DUR-A3-2A-A	DUR-A3-2A-A	DURA-CUF, LARGE ADULT, DINAClick, 31 - 40 CM, 80369-5, 5/ BOX Specific Labeling

Table 13 SOFT-CUF® blood pressure cuffs, dual-tube, GEMS-IT

GE HealthCare Part Number	Manufacturer Part Number	NIBP Accessory Description
SFT-N1-2B	SFT-N1-2B	SOFT-CUF, Cuff, Neonatal type 1, 2T Neo-Snap
SFT-N2-2B	SFT-N2-2B	SOFT-CUF, Cuff, Neonatal type 2, 2T Neo-Snap
SFT-N3-2B	SFT-N3-2B	SOFT-CUF, Cuff, Neonatal type 3, 2T Neo-Snap
SFT-N4-2B	SFT-N4-2B	SOFT-CUF, Cuff, Neonatal type 4, 2T Neo-Snap
SFT-N5-2B	SFT-N5-2B	SOFT-CUF, Cuff, Neonatal type 5, 2T Neo-Snap
SFT-P1-2A-5S	SFT-P1-2A-5S	SOFT-CUF, INFANT, DINACLICK, 08 - 13 CM, 80369-5, 20/ BOX
SFT-P2-2A-5S	SFT-P2-2A-5S	SOFT-CUF, CHILD, DINACLICK, 12 - 19 CM, 80369-5, 20/ BOX
SFT-P2-2A-L-5S	SFT-P2-2A-L-5S	SOFT-CUF, CHILD LONG, DINACLICK, 12 - 19 CM, 80369-5, 20/ BOX
SFT-P2-2A-INT-5S	SFT-P2-2A-INT-5S	SOFT-CUF, CHILD, DINACLICK, 12 - 19 CM, INTERNATIONAL, 80369-5, 20/ BOX
SFT-A1-2A-5S	SFT-A1-2A-5S	SOFT-CUF, SMALL ADULT, DINACLICK, 17 - 25 CM, 80369-5, 20/ BOX
SFT-A1-2A-L-5S	SFT-A1-2A-L-5S	SOFT-CUF, SMALL ADULT LONG, DINACLICK, 17 - 25 CM, 80369-5, 20/ BOX
SFT-A1-2A-INT-5S	SFT-A1-2A-INT-5S	SOFT-CUF, SMALL ADULT, DINACLICK, 17 - 25 CM, INTERNATIONAL, 80369-5, 20/ BOX
SFT-A2-2A-5S	SFT-A2-2A-5S	SOFT-CUF, ADULT, DINACLICK, 23 - 33 CM, 80369-5, 20/ BOX
SFT-A2-2A-INT-5S	SFT-A2-2A-INT-5S	SOFT-CUF, ADULT, DINACLICK, 23 - 33 CM, INTERNATIONAL, 80369-5, 20/ BOX
SFT-A2-2A-L-5S	SFT-A2-2A-L-5S	SOFT-CUF, ADULT LONG, DINACLICK, 23 - 33 CM, ENG., 80369-5, 20/ BOX
SFT-A2-2A-L-INT-5S	SFT-A2-2A-L-INT-5S	SOFT-CUF, ADULT LONG, DINACLICK, 23 - 33 CM, INTERNATIONAL, 80369-5, 20/ BOX
SFT-A3-2A-5S	SFT-A3-2A-5S	SOFT-CUF, LARGE ADULT, DINACLICK, 31 - 40 CM, 80369-5, 20/ BOX
SFT-A3-2A-INT-5S	SFT-A3-2A-INT-5S	SOFT-CUF, LARGE ADULT, DINACLICK, 31 - 40 CM, INTERNATIONAL, 80369-5, 20/ BOX
SFT-A3-2A-L-5S	SFT-A3-2A-L-5S	SOFT-CUF, LARGE ADULT LONG, DINACLICK, 31 - 40 CM, 80369-5, 20/ BOX

Table 14 CLASSIC-CUF® blood pressure cuffs, dual-tube, GEMS-IT

GE HealthCare Part Number	Manufacturer Part Number	NIBP Accessory Description
CLA-N1-2B	CLA-N1-2B	CLASSIC-CUF, Cuff, Neonatal type 1, 2T Neo-Snap
CLA-N2-2B	CLA-N2-2B	CLASSIC-CUF, Cuff, Neonatal type 2, 2T Neo-Snap
CLA-N3-2B	CLA-N3-2B	CLASSIC-CUF, Cuff, Neonatal type 3, 2T Neo-Snap
CLA-N4-2B	CLA-N4-2B	CLASSIC-CUF, Cuff, Neonatal type 4, 2T Neo-Snap
CLA-N5-2B	CLA-N5-2B	CLASSIC-CUF, Cuff, Neonatal type 5, 2T Neo-Snap
CLA-P1-2A-5S	CLA-P1-2A-5S	CLASSIC-CUF, INFANT, DINACLICK, 08 - 13 CM, 80369-5, 20/ BOX
CLA-P2-2A-5S	CLA-P2-2A-5S	CLASSIC-CUF, CHILD, DINACLICK, 08 - 13 CM, 80369-5, 20/ BOX
CLA-A1-2A-5S	CLA-A1-2A-5S	CLASSIC-CUF, SMALL ADULT, DINACLICK, 17 - 25 CM, 80369-5, 20/ BOX

Table 14 CLASSIC-CUF® blood pressure cuffs, dual-tube, GEMS-IT (Table continued)

GE HealthCare Part Number	Manufacturer Part Number	NIBP Accessory Description
CLA-A2-2A-5S	CLA-A2-2A-5S	CLASSIC-CUF, ADULT, DINACLICK, 23 - 33 CM, 80369-5, 20/ BOX
CLA-A2-2A-L-5S	CLA-A2-2A-L-5S	CLASSIC-CUF, ADULT LONG, DINACLICK, 23 - 33 CM, 80369-5, 20/ BOX
CLA-A3-2A-5S	CLA-A3-2A-5S	CLASSIC-CUF, LARGE ADULT, DINACLICK, 31 - 40 CM, 80369-5, 20/ BOX
CLA-A3-2A-L-5S	CLA-A3-2A-L-5S	CLASSIC-CUF, LARGE ADULT LONG, DINACLICK, 31 - 40 CM, 80369-5, 20/ BOX
CLA-P1-2A-X-5S	CLA-P1-2A-X-5S	CLASSIC-CUF ISO, INFANT, DINACLICK, 08 - 13 CM, 80369-5, 20/ BOX
CLA-P2-2A-X-5S	CLA-P2-2A-X-5S	CLASSIC-CUF ISO, CHILD, DINACLICK, 12 - 19 CM, 80369-5, 20/ BOX
CLA-A1-2A-X-5S	CLA-A1-2A-X-5S	CLASSIC-CUF ISO, SMALL ADULT, DINACLICK, 17 - 25 CM, 80369-5, 20/ BOX
CLA-A2-2A-X-5S	CLA-A2-2A-X-5S	CLASSIC-CUF ISO, ADULT, DINACLICK, 23 - 33 CM, 80369-5, 20/ BOX
CLA-A3-2A-X-5S	CLA-A3-2A-X-5S	CLASSIC-CUF ISO, LARGE ADULT, DINACLICK, 31 - 40 CM, 80369-5, 20/ BOX

Table 15 ONE-CUF® Blood Pressure Cuffs, Dual-Tube, DinaClick Connector

GE HealthCare Part Number	Manufacturer Part Number	NIBP Accessory Description
ONE-P1-2A	ONE-P1-2A	ONE-CUF, INFANT, DINACLICK 80369-5, 08 - 13 CM, 20/BOX
ONE-P2-2A	ONE-P2-2A	ONE-CUF, CHILD, DINACLICK 80369-5, 12 - 19 CM, 20/BOX
ONE-A1-2A	ONE-A1-2A	ONE-CUF, SMALL ADULT, DINACLICK 80369-5, 17 - 25 CM, 20/BOX
ONE-A2-2A	ONE-A2-2A	ONE-CUF, ADULT, DINACLICK 80369-5, 23 - 33 CM, 20/BOX
ONE-A2-2A-L	ONE-A2-2A-L	ONE-CUF, ADULT LONG, DINACLICK 80369-5, 23 - 33 CM, 20/BOX
ONE-A3-2A	ONE-A3-2A	ONE-CUF, LARGE ADULT, DINACLICK 80369-5, 31 - 40 CM, 20/BOX

Temperature accessories

Table 16 Welch Allyn accessories

GE HealthCare Part Number	Manufacturer Part Number	Temperature Accessory Description
2068729-002	02895-000	WELCH ALLYN ORAL PROBE - 9FT CORD WITH PROBE WELL
2068729-003	02895-100	WELCH ALLYN RECTAL PROBE - 9FT CORD WITH PROBE WELL
2068729-004	05031-101	WELCH ALLYN PROBE COVERS, DISPOSABLE, 1000 COUNT, 25 PER BOX
2068729-005	05031-110	WELCH ALLYN PROBE COVERS, DISPOSABLE, 10000 COUNT, 25 PER BOX
2068729-006	06138-000	WELCH ALLYN CALIBRATION KEY
2068729-007	01802-110	WELCH ALLYN 9600 PLUS CALIBRATION TESTER

Table 17 Exergen accessories

GE HealthCare Part Number	Manufacturer Part Number	Temperature Accessory Description
EX134202	EX134202	134202 CAP DISPENSER (EXERGEN)
EX134203	EX134203	RESPOSABLE CAPS, EXERGEN, BOX/1000
EX134205	EX134205	RESPOSABLE CAPS, EXERGEN, BOX/5000
EX129462	EX129462	Disposable tubular protective sheath, 250/box
EX129003	EX129003	129003 CALIBRATION VERIFCTN KIT, EXERGEN
EX129012	EX129012	Heat Generator, Exergen
5921582	124227-AF-GE	Sensor, Exergen TAT, Arterial REF, °F, USB, white (TAT5000S-USB-GE2)
5921583	124227-AC-GE	Sensor, Exergen TAT, Arterial REF, °C, USB, white ((TAT5000S-USB-GE2)
5921584	124227-OF-GE	Sensor, Exergen TAT, Oral REF, °F, USB, white (TAT5000S-USB-GE2)
5921585	124227-OC-GE	Sensor, Exergen TAT, Oral REF, °C, USB, white (TAT5000S-USB-GE2)
2044860-013	124260	Sensor, Exergen TAT, Arterial REF, °F, USB, blue (TAT5000S-USB-GE)
2044860-014	124261	Sensor, Exergen TAT, Arterial REF, °C, USB, blue (TAT5000S-USB-GE)
2044860-015	124262	Sensor, Exergen TAT, Oral REF, °F, USB, blue (TAT5000S-USB-GE)
2044860-016	124263	Sensor, Exergen TAT, Oral REF, °C, USB, blue (TAT5000S-USB-GE)

Table 18 HeTaida accessories

GE HealthCare Part Number	Manufacturer Part Number	Temperature Accessory Description
5919859	HTD8808C (CE)	HeTaida Non Contact Infrared Body Thermometer (CE)
5919860	HTD8808C (FDA)	HeTaida Non Contact Infrared Body Thermometer (FDA)
5919862	HTD8808C (CN)	HeTaida Non Contact Infrared Body Thermometer (China)

Recorder paper

Table 19 Recorder paper, GRAPHIC CONTROLS

GE HealthCare Part Number	Manufacturer Part Number	Accessory Description
2106823-001	32028978	Thermal paper 50 mm x 27.4 M, Blank, 48 rolls

Supplies

Mountings

- GE HealthCare devices are designed to allow for reliable mounting attachment to a mount. The mounting attachment is approved for use with the listed mounts.
- Follow mount manufacturer instructions for installation and loading.

Table 20 GCX mounting

GE HealthCare Part Number	Manufacturer Part Number	Accessory Description
2033261-001	QCT-0001-30	COUNTERTOP MOUNT BASE 3.2"HEIGH 4"X4"
2002720-001	HZ-0001-01	Low Profile Horizontal Channel Mount w/ Slide-In Mounting Plate
2002746-002	ME-0047-13	GEHC 15" Flat Panel VHM Arm - No Channel
407349-009	WMM-0002-02	M Series 12"/30.5 cm Pivot Arm with Slide-In Mounting Plate
407349-022	WMM-0002-06	M Series 12 x 12"Articulating Arm with Slide-In Mounting Plate
2038478-001	WMM-0006-04D	M SERIES ARM UTILITY HOOK
416537-006	WC-0002-07	Seismic PolyMount 25"/63.5 cm Channel
416537-012	WC-0002-09	Seismic PolyMount 31"/79 cm Channel
411959-001	WC-0002-05	CHANNEL, SEISMIC 19 INCH GCX CORP
5857340	DX-0044-16	Bed Mounting for 10 inch_12 inch monitor

Table 21 Colytech mounting

GE HealthCare Part Number	Manufacturer Part Number	Accessory Description
5838484	AMF0001	FLA SERIES PIVOT ARM 12 INCH
5838485	AY.00005.PA	FLA SERIES UTILITY CABLE HOOK
5838500-02	AY.00002.AS	MOUNT, WALL CHANNEL 19 INCH
5838500-03	AY.00003.AS	MOUNT, WALL CHANNEL 25 INCH
5838500-04	AY.00004.AS	MOUNT, WALL CHANNEL 31 INCH
5937979	5937979	WALL MOUNTING
6660123-60	6660123-60	5 Inch Mounting plate

Other supplies

Table 22 Barcode scanner

GE HealthCare Part Number	Manufacturer Part Number	Accessory Description
5813075	CR1500-K101-C508F-C241	COMPACT USB BARCODE SCANNER AND 8FT/2.4M COILED CABLE, CODE
5942351	1470G2D-5USB-MS	Honeywell 1470G USB Barcode scanner, 3m coiled cable

Table 23 Roll stand

GE HealthCare Part Number	Manufacturer Part Number	Accessory Description
5923904	CTF0011	Roll Stand
6660123-71	GE0682AS1Y3A	Roll Stand bottom basket Kit (1/box)
6660123-72	GE0683AS1Y3A	Roll Stand bottle holder Kit (1/box)
6660123-73	GE0680AS1Y3A	Roll Stand holder Kit (75mm height, 5/box)

Table 23 Roll stand (Table continued)

GE HealthCare Part Number	Manufacturer Part Number	Accessory Description
6660123-74	GE0681AS1Y3A	Roll Stand holder Kit (50mm height, 5/box)
6660123-75	GE0686AS1Y3A	Roll Stand cable holder kit (5/box)

Table 24 Active RFID

GE HealthCare Part Number	Manufacturer Part Number	Accessory Description
2105601-001	GE-MB1000-25-WR	GE ASSET MOBILE BEACON - 25 PACK