

PRODUCT DATASHEET TROLLEY 2.101126



SCHILLER
The Art of Diagnostics

Overall dimensions and weight:	515 mm (W) × 515 mm (D) × 1125 mm (H) - 14.3 kg
Top Console in steel with handle in abs	
Work Height	1055 mm (center of the device)
Dimensions of the device support	242 mm (W) × 167 mm (D)
Dimensions with handle	326 mm (W) × 360 mm (D)
Aluminium Column with support	
Base footprint	515 × 515 mm
4 wheels with brake	ø 100 mm
Power supply support in steel	
Dimensions	116 mm (W) × 184 mm (D) × 39 mm (H)
Inner Dimensions	60 mm (W) × 180 mm (D) × 37 mm (H)
Basket in Polycarbonate	
Dimensions	350 mm (W) × 250 mm (D) × 105 mm (H)
Height	Adjustable
Load limit	5 Kg
3 Pass-throughs with grommets in abs	for cable management
Overall load limit	30 Kg
Accessories	Holder barcode scanner
	ECG Cable support
	14 -wire patient cable support
	Basket
	Drawer
	Holder SpiroScout SP plus
	DT-80 holder
	DT-100 holder
Packaging dimensions and weights	Basis 2.101129: 800 mm (W) × 600 mm (D) × 300 mm (H) - 12.7 kg Console 2.101133: 350 mm (W) × 340 mm (D) × 115 mm (H) - 3.6 kg
CE Declaration of Conformity	EN 60601-1
	Directive 93/42/CE medical device



SCHILLER
A M E R I C A S

Call us: **1-888-845-8455**
Visit us at **schillerservice.com**
10903 NW 33rd St, Doral, FL 33172, USA