

SILVERFOX

Professional design, comfortable experience

Stool

8C01



Silverfox Corporation Limited

Tel : 0086-750-2633320 2633805
sales@silverfox.cn www.silverfoxgp.com

Address : No.18,1st TongLe Road,TangXia Town,PengJiang
District,JiangMen City,GuangDong Province,China



Relax mode

The Relax mode keep the body in a relaxed state. Lift the backrest to the level of dorsal and tilt backward, adjust the seat-rest at horizontal angle.



Consulting mode

The Consulting mode facilitates face to face communication between practitioner and patient in a comfortable posture.

Adjust the backrest to upright position, and seat-rest to horizontal angle, this posture helps support the spine at a straight and comfortable state.



Operation mode

The Operation mode is specially suitable for working at a forward leaning position.

Adjust the backrest to support the waist, and the seat-rest at incline angle. Helping the thigh and calf to form a 120-degree angle will reduce stress of the spine and prevent injury.

Specifications

Net weight: 10.7 kg / 23lbs

Product size : 600x600x750mm

Gross weight: 12.6 kg / 27.2lbs

Package Size : 590x300x600mm



470-590mm

600mm



Medical



Simple